


FUJI MEDICAL DRY IMAGING FILM
DI-AT


Fuji Medical Dry Imaging Film DI-AT is designed for use with FUJI MEDICAL DRY IMAGER DRYPIX 1000/3000.

■ **SENSITIVITY & CONTRAST**

Sensitivity and contrast of the DI-AT are suitably designed for thermal-recording with FUJI MEDICAL DRY IMAGER DRYPIX 1000/3000.




■ **FILM STRUCTURE**



■ **AVAILABLE SIZE & QUANTITY**

Film Size	Quantity / Package
① 25.7 × 36.4 cm	100 sheets + 1 cleaning sheet
② 35.4 × 43.0 cm	100 sheets + 1 cleaning sheet

■ DIMENSIONS**■ STORAGE CAUTIONS**

- Make sure to keep unrecorded film in its moisture-proof bag, stored in a cool, dark place (23°C/74°F or below).
- Observe the following points in handling and storing recorded film.
Store recorded film in a cool, dry place of low temperature and low humidity. The higher the temperature and humidity, the more the density of recorded images will increase. Long-term storage at high temperatures, high humidity and/ or daylight room conditions may cause discoloration.

Results of accelerated testing indicate the following life spans for recorded film, in terms of how long it takes for a 10% rise in density to occur in the places in which image density is 1.1 at time of output.

① When stored at 25°C, 60%RH, approx. 30 years.

② When stored at 30°C, 60%RH, approx. 5 years.

If storage temperature exceeds these levels, lifespan will be as follows.

③ When stored at 35°C, 60%RH, approx. 6 months.

④ When stored at 40°C, 60%RH, approx. 1 week.

Thus, make sure not to leave recorded film in places exposed to direct sunlight, next to heating equipment, or inside automobiles.

Under the following environmental conditions, the lifespan will be approx. 5 years.

Whithin a one-year period:

- When stored 4 months at 35°C, 60%RH for 8 hours a day and 28°C, 60%RH for 16 hours a day.
- When stored 8 months at 28°C, 60%RH, 24 hours a day.

*For other details on handling of DI-AT film, refer to the operation manual.