

Unsurpassed Performance in Medical Radiography

The Choice for Versatility, Reliability and Value



AM3 Integrated Table System

- **Digital imaging compatible**

AMERICOMP

Affordable Choices for General Radiography.

Since 1987, the quality standards of Americomp™ have been unsurpassed. Americomp™ systems remain the affordable choice in today's marketplace, offering the same long-term dependable performance that has been synonymous with the Americomp™ name. Our five-year parts warranty secures your investment today and guarantees your investment value tomorrow.



AM3 Integrated Table System Features:

- Tubestand integration to table saves valuable space and does not require wall or ceiling modifications
- Unique off-table option provides a full range of off-table procedures including standing knee and ankle exams (shown above).
- Versatile tube positioning offers a wide range of procedures
- Flat, extra-wide, table top allows easy patient positioning
- Heavy-duty construction with 300 lb. patient rating
- Wallstand provides effortless positioning of stationary grid cabinet or optional reciprocating bucky

High Performance Generators

Complement any Americomp radiographic system with a high-quality, state-of-the-art generator. Contact your Americomp dealer for details.



AM4 Stationary Table System Features:

- Floor-to-wall or ceiling tubestand offers versatile tube positioning including off-table standing knee and ankle exams
- Flat, extra-wide, table top allows easy patient positioning
- Heavy-duty construction with 300 lb. patient rating
- Wallstand provides effortless positioning of stationary grid cabinet or optional reciprocating bucky

Digital-ready Systems

Americomp radiographic equipment is completely compatible with most digital imaging systems. Digital imaging technology may be added at any time.



AM5 Compact Table System Features:

- Perfect for "no space" situations
- Fits into rooms measuring as small as 7 ft. x 9 ft.
- Floor-to-wall or ceiling tubestand offers versatile tube positioning
- Flat table top allows easy patient positioning
- Heavy-duty construction with 300 lb. patient rating
- Wall mounted cassette holder saves space and does not require handling of cassette tray

#QEC15774
ISO 9001:2000 Certified

