

QT-710

Radiographic Mobile Fixed Table

650 lb. Patient
Weight Capacity



QT-710 MOBILE FIXED TABLE

Quantum's QT-710 Mobile Fixed Table offers mobile flexibility, fast positioning and ease-of-use. This unique ability makes the QT-710 Table an excellent choice for simplicity in positioning patients for exams, and is excellent for performing supine exams. The table boasts large 5" locking castor wheels, which have the ability to be locked individually depending on use. The locking feature maintains table stability and steadiness in any imaging environment. The large castor wheels ensure smooth table movement with a patient weight capacity up to 650 lbs (272 kg). The completely flat table-top design is low maintenance and easy to clean. Made of a robust low absorption radiolucent material, the table-top is durable.

Specifications

- Patient Weight Capacity: 650 lbs. (272 kg)
- Table Length: 74" (188 cm)
- Table Width: 28" (71 cm)
- Table Height: 31.75" (81 cm)
- Table Length: 74" (188 cm)
- Includes one set (2) Patient Grips
- Low absorption tabletop material
- Locking 5" castor wheels, individually locking

Specifications subject to change without prior notice.

