

# 3308B

## ROTATING TUBE STAND




### Technical Specifications:

Ceiling height:	Minimum 96in (244cm) Maximum 116in (295cm)
Vertical travel:	Minimum 12in (30.5cm) Maximum 75in (190.5cm)
Longitudinal travel:	72in (183cm) with 8ft rail
Longitudinal rails length:	8ft standard
Transverse travel:	9in (23cm)
Column rotation:	±180° with detents 0°, 90°, 90°
Tube rotation:	+180°, -180°
ETL Listed	

### Tube Stand Options:

- 140KHU, 200KHU or 300KHU X-ray tubes
- Manual, semi-automatic or fully automatic collimators
- Trunnion mount in place of standard platform mount
- 12ft, 10ft or 6ft longitudinal rails



**DEL MEDICAL**

Toll Free: (800) 800-6006  
 Fax: (847) 288-7011  
 Email: info@delmedical.com  
 www.delmedical.com  
 241 Covington Drive, Bloomingdale, IL 60108

**Universal**

Universal branded products are manufactured and distributed by UMG/Del Medical. Del Medical is an ISO 13485 & ISO 9001 Certified facility. Specifications are subject to change without notice.

**UMG**

Toll Free: (800) 261-9808  
 Fax: (914) 835-6111  
 Email: info@umgxray.com  
 www.umgxray.com  
 28 Calvert Street, Harrison, NY 10528

08.13 #9000-3308BUX

# 3308B

## ROTATING TUBE STAND

The Universal 3308B rotating tube stand offers the versatility needed to facilitate a wide range of general purpose radiographic procedures. The sturdy tube mount assembly ensures stability and precision alignment.



Longitudinal lock, roll lock and vertical lock release button with positive indexing at 90 degrees for fast, accurate cross table positioning.

### Key Features

- Pressure activated column rotation for fast, accurate positioning
- Electric vertical and longitudinal locks provide maximum holding power
- 180° column rotation and centering detent allow for rapid patient positioning
- 9in (23cm) of transverse travel for total table top accessibility

