

VERTI-Q Wall Stands

**All of the VERTI-Q Wall Stands can also be configured with a choice of multi-speed buckies, fixed grids or automatic cassette loading systems for use with traditional film systems.*

QW-420 "VERTI-Q" Wall Stand

- Features the exclusive "EZ-Glide" Hand control for easy and precise movement, grip rotates +105°
- Custom enclosed frame for attractive appearance includes patient chin rest
- Low absorption front cover material with cassette and AEC Indicators
- Vertical Travel: 60.0 (±1.0)" (1524.0 mm) with a 15.0 (±0.5)" (381.0 mm) minimum Focal Spot-to-Floor Distance
- FAIL-SAFE electromagnetic braking system plus integral counterbalancing, ensure safe and easy use
- Accepts cassette sizes of up to 14" x 17" (35 x 43 cm)

NOTE: Specify either Right or Left Hand Loading

**Specifications are subject to change*



QW-420-S "VERTI-Q SIDE MOUNTED" Space Saving Wall Stand

Same as QW-420, except image receptor is mounted to the side of the structure

- Features the exclusive "EZ-Glide" Hand control for easy and precise movement, grip rotates +105°
- Custom enclosed frame for attractive appearance includes patient chin rest
- Low absorption front cover material with cassette and AEC Indicators
- Vertical Travel: 60.0 (±1.0)" (1524.0 mm) with a 15.0 (±0.5)" (381.0 mm) minimum Focal Spot-to-Floor Distance
- FAIL-SAFE electromagnetic braking system plus integral counterbalancing, ensure safe and easy use
- Accepts cassette sizes of up to 14" x 17" (35 x 43 cm)

NOTE: Specify either Right or Left Hand Loading

NOTE²: QW-HG30 Overhead handgrip is available on Right Loading Only

NOTE³: QW-HG20 Side Mounted Handgrips & QW-EX4.5 Extension (4.5") are *NOT* available on the QW-420-S

**Specifications are subject to change*



QW-420-MT "VERTI-Q TILT" Manual Tilting Wall Stand

- Features the exclusive "EZ-Glide" Hand control for easy and precise movement, grip rotates +105°
- Custom enclosed frame for attractive appearance includes patient chin rest
- Low absorption front cover material with positioning and AEC Indicators
- Vertical Travel: 57.0 (±1.0)" (1447.8 mm) with a 15.0 (±0.5)" (381.0 mm) minimum Focal Spot-to-Floor Distance
- Manual Tilting Movement (+90°/-15°), with stops at -15° +90° detents at 0°
- FAIL-SAFE electromagnetic braking system plus integral counterbalancing, ensure safe and easy use
- Accepts cassette sizes of up to 14" x 17" (35 x 43 cm)

NOTE: Specify either Right or Left Hand Loading

**Specifications are subject to change*



QW-420-T "VERTI-Q TILT" Motorized Tilting Wall Stand

- Features the exclusive "EZ-Glide" Hand control for easy and precise movement, grip rotates +105°
- Custom enclosed frame for attractive appearance includes patient chin rest
- Low absorption front cover material with positioning and AEC Indicators
- Vertical Travel: 60.0 (±1.0)" (1524.0 mm) with a 15.0 (±0.5)" (381.0 mm) minimum Focal Spot-to-Floor Distance
- Motorized Tilting Movement (+90 / -15°), with stops at 0 and 90 degrees
- FAIL-SAFE electromagnetic braking system plus integral counterbalancing, ensure safe and easy use
- Accepts cassette sizes of up to 14" x 17" (35 x 43 cm)

NOTE: Specify either Right or Left Hand Loading

NOTE²: If Power Supply is Required - PS-8 (AY20-084) Must Be Ordered Separately

**Specifications are subject to change*



OPTIONS FOR "VERTI-Q" WALL STANDS

QW-HG20 Patient "Side Mounted" Handgrips: (One Set), Mounts To "VERTI-Q" Wall Stand

QW-HG30 Patient "Overhead" Handgrip: Mounts To "VERTI-Q" Wall Stand

- Grip can be rotated, with lock positions at 0°, 90° and 180° (grip may be removed, for certain positions)

QW-EX8 Extension Bracket (8"): Extends Bucky front cover an additional 8" from Tubestand

- NOTE: *ONLY AVAILABLE ON "QW-420 & QW-420-D"* (aids in positioning for wheelchair patients)

QW-EX4.5 Extension (4.5"): Extends Bucky front cover an additional 4.5" from Tubestand

- NOTE: *ONLY AVAILABLE ON "QW-420-MT, QW-420-MTD, QW-420-T & QW-420-TD"* (for alignment of various ceiling structures, to enable "center line to x-ray beam")