

CMP 200®DR series

Radiographic X-Ray generators up to 400 kHz – unrivalled reliability and performance.

4 different models to choose from:

- 40 kW
 - 50 kW
 - 65 kW
 - 80 kW
- Digital Radiography Interfaces*
 - Dual Speed Starter*
 - Cost-effective for all radiographic applications (hospitals, medical clinics, chiropractic and veterinary applications)
 - Field upgradable for your future imaging requirements.

Small. Compact. Lightweight. Cost Effective.

* Consult factory for details and other options

X-RAY
Generators of Choice

CPI family of 'Generators-of-



Membrane Console Standard

'Consoles-of-

ly of X-Ray rs-of-Choice'.



Touchscreen Optional

es-of-Choice'

Look at these Features

- CMP200[®] DR generator utilizes CPI's proven up to 400 kHz⁺ inverter technology and software architecture.
- 80 kW enhanced features and performance to meet the most demanding requirements of today's busy imaging centers.
- Anatomical programs available in many languages that are easily edited via the Operator Control Console
- Technique Selection:
 1. kV with AEC*
 2. kV/mAs
 3. kV/mA/Time
 4. Patient/Body Part Thickness
(Membrane Control Console)
- Low kV ripple equal to constant potential performance maximizes image quality while reducing dose.
- Automatic Tube Calibration
- GenWare[®] Service Software allows fast and easy set up as well as system diagnostics for off-site analysis
- Collimator Power Supply
- Basic Tomography Interface

Options:

- Digital Radiography Interfaces*
- Touchscreen Control Console (with embedded GenWare[®] Service Software)
- Dual Speed Starter for higher power
- Exposure Hand Switch
- Mini Console (for Integrated Systems)
- Console Wall Mount
- Automatic Exposure Control (AEC) Interface
- Dose Area Product (DAP) Interface
- Auto Positioning Interface

Compliance with FDA, UL/CSA, MDD, CE, CQC, and SFDA Requirements

* Consult factory for details and other options



Exposure Hand Switch Optional



Mini Console Optional



Wall Mount
Optional



Membrane Console Pedestal
Optional

World Class Product Support

One of our greatest strengths is our world class reputation for customer support. A dedicated support team is only a phone call away at 905-877-0161 – or at our North American toll-free number 1-888-CPI-XRAY (274-9729)

Call us. We are here to help you with all of your installation needs.



Communications & Power Industries Canada Inc.
communications & medical products division
45 River Drive / Georgetown, Ontario / Canada L7G 2J4
Hot Line Telephone: 1-888-CPI-XRAY
TEL: 905-877-0161 / FAX: 905-877-5327
E-MAIL: marketing@cmp.cpii.com
www.cpii.com



Some features shown may be optional. Please refer to Data Sheet.