

CR Pro

Advanced CR for Accuracy and Productivity

Overview

The CR Pro is a complete and fully automated solution that addresses all of your digital x-ray needs. It integrates seamlessly with your existing practice and has the lowest cost of ownership and best customer support in the industry.

Made in USA

All Radlink products are designed and manufactured in Southern California at our high tech manufacturing plant. Our technical support team is also all located in the US to respond to any issues you may have.

Complete System

The CR Pro is a compact mobile unit with embedded CPU and complete image acquisition and DICOM transmission capabilities. Image Acquire, DICOM Print, Multi-Modality Router, Modality Worklist, Local Archive, CD burner, and touch-screen flat panel display are all included.

User Friendly

Our touch-screen software interface is intuitive and user-friendly. Plate processing is rapid (less than 45 seconds) for improved productivity. Installation and training can be completed in a single day and will fit into any workflow environment - whether stand-alone, PACS, or teleradiology. The CR Pro is a self-contained unit on wheels that can be easily moved throughout your facility.

Superior Image Quality

CR Pro produces image quality that is top of the line. Our core technology utilizes proprietary sealed fibers (patent pending) that maximize energy transfer and preserve image quality. Radlink's advanced post-processing software, included at no extra charge, generates images with improved contrast and visibility even at lower x-ray doses.

Lowest Cost of Ownership

Our cassette-based system maximizes the life of costly phosphor-image plates. The operator never risks causing damage by handling the image plate. To further maximize service life, the image plate is mounted on a rigid substrate and is never bent or touched by the feed mechanism. On top of this, a typical private practice, with 20-30 x-rays per day, saves about \$20,000 per year on film costs. The CR Pro pays for itself in less than 2 years!

Improves Patient Care

With the CR Pro your images are ready to view within seconds, leading to quicker diagnosis and reduced waiting time for patients. The web-based PACS system enables remote diagnosis through teleradiology, and the viewing software helps illustrate and explain complex diagnoses to the patient.

Unparalleled Customer Support

Every CR Pro is backed by our one-year parts and service warranty, which can be extended up to five years. Unlike our competitors who require you to go through multiple emails, we provide immediate direct phone access to our knowledgeable technical support staff. If we can't solve the problem remotely we'll send someone to your site the next business day so you can get back to the business of medicine.

