




KONICA MINOLTA

ImagePilot M – multi-user, feature rich all-in-one

The All-In-One digital radiography solution configured for most private practice healthcare providers who want the simplest transition to digital and room for growth.



Standard ImagePilot Features:

- All-In-One Design—Easiest for the staff to learn. Reduces IT support requirements.
- AutoPilot—Automatic image processing for unsurpassed consistency in image quality.
- Simplified Image Acquisition—Single click workflow increases productivity.
- Same User Interface At Every Workstation—learn only one system.
- Easy-To-Use And Feature Rich Viewing Screen—Image adjustments, annotations, measurements, inputs and outputs, customizable for basic to advanced users.

ImagePilot M highlights:

- Excellent Speed—60 cassettes per hour upgradable to 76.
- Redundant Storage—mirrored 500GB hard drives (RAID 1).
- Diagnostic Quality Display—Review images with confidence on a calibrated diagnostic monitor.
- Distribute Workflow with Multiple Users—Three additional licenses (4 total inc. the server).
- Import Priors or Export Patient Images to CD—Better than film for transportability and viewing.
- Power Protection—UPS for the server to provide power during outages or protect against anomalies.

ImagePilot

The essentials of imaging

ImagePilot Features and Specifications

Hardware Components	ImagePilot LE	ImagePilot LS	ImagePilot M	ImagePilot PRO
Imaging Station				
Server Platform	Desktop	Tower	Tower	Tower
Hard Drive Capacity	250GB (1x250GB)	500GB (2x500GB RAID 1)	500GB (2x500GB RAID 1)	1.5TB (3x750GB RAID 1)
RAM	2GB	2GB	2GB	4GB
Display Options	19" Lenovo 21" or 24" Eizo	19" Lenovo 21" or 24" Eizo	21" or 24" Eizo	21" or 24" Eizo
UPS for Imaging Station	Optional	Optional	√	√
REGIUS 110 Series				
CR Reader	REGIUS 110	REGIUS 110	REGIUS 110	REGIUS 110
Throughput (14x17 cassettes per hour)	60 (upgradeable to 76)	40	60 (upgradeable to 76)	60 (upgradeable to 76)
Backup/Extended Storage				
DVD Backup Tool- For Disaster Recovery	√	√	√	√
Backup NAS, 2TB or 4TB RAID 5	Optional	Optional	Optional	Included, 2TB RAID 5
Extended Storage NAS, 2TB or 4TB RAID 5	Optional	Optional	Optional	Optional
Software Features				
Patient Registration				
Manual - data entry	√	√	√	√
Automatic-Worklist Toolkit (includes DICOM MWL)	√	Optional	Optional	√
Image Input				
AutoPilot Image Processing for CR Image Acquisition	√	√	√	√
Import images from CD/DVD, USB/folder-DICOM/JPEG	Optional	Optional	√	√
Additional DICOM Modality Acquisition	Optional	Optional	Optional	Optional
Image Review Tools				
View, Annotate, Measure	√	√	√	√
Report Tool (requires office)	Optional	Optional	√	√
Simple Note Tool	√	√	√	√
Image Output				
Export images to CD/DVD w/ DICOM viewer or to folder	√	√	√	√
Multi-patient Export to CD/DVD or files	Optional	Optional	Optional	Optional
Print - Windows Printing	√	√	√	√
DICOM Store PACS (first target)	√	√	√	√
Distributed Workflow				
Concurrent Users (including server)	2	3	4	5
Max # of authorized client computers	19	19	19	19
Remote Viewing (IP client with compressed previews)	Optional	Optional	Optional	Optional
Advanced Connectivity				
Worklist Toolkit - DICOM MWL	√	Optional	Optional	√
DICOM Query Retrieve (SCP or SCU)	Optional	Optional	Optional	Optional
EMR Integration (i.e. HL7)	Optional	Optional	Optional	Optional

Features and Specifications are subject to change without notice.

