


TRANSFER SETS



TRANSFER SETS:


Sterile transfer sets with various designs to fill syringes from multiple sizes and styles of containers.

TRANSFER SETS

CONTRAST SPIKE IS AVAILABLE WITH BALL CHECK TO PREVENT LEAKAGE



LONG AND SHORT SPIKES AVAILABLE



SUPPORTS FASTER FILL RATES

AVAILABLE WITH SWABABLE VALVE



- 1 lower cost
- 2 enhanced productivity
- 3 supports optimum clinical outcomes



www.coeurinc.com

1.800.296.5893
1.615.547.7923
Coeur, Inc. Corporate
704 Cadet Court
Lebanon, TN 37087