

# Medlink CR-X5 System

High-End Versatile CR System which offers an ideal solution for any private practice or decentralized CR environment.

Medlink CR-X5



# MEDLINK CR-X5 SYSTEM



MEDLINK IMAGING'S COMPACT CR-X5 DIGITIZER IS AN IDEAL SOLUTION FOR SMALL PRIVATE PRACTICES, OFFERING A COST-EFFECTIVE SOLUTION INTO DIGITAL RADIOGRAPHY WITHOUT COMPROMISING ON IMAGE QUALITY AND PATIENT WORKFLOW.

The CR-X5 has a low cost of ownership thanks to its modular design and ease of installation; as well as, being very easy to use and maintain. The CR-X5 Solution also includes Datalink, a Fully Managed HIPAA Compliant Disaster Recovery Service; along with, an Integrated Mini-PACS Offering with the ability to view and analyze images from any location thanks to the system's local and web viewing licenses.

## CR-X5 CR System 14x17 Plates Per Hour

40-pph @ 100µm / 50-pph @150µm / 60-pph 200µm

Multiple Sized Cassette Feed (14x17, 24x30, 18x24 and 5x30)  
Included are Two (2) MD 1.0 14x17 (35x43) Cassettes

Data Acquisition: 20 bits/pixel



Resolution: 3480x4248 (14.8MP)

Output to Processor: 16 bits/pixel

Power: 120-240 VAC, 60 Hz

Weight: 63 lbs.

Dimensions: 22.83 w x 27.56 d x 18.54 h



## CR-X5 Acquisition Workstation

1TB in a RAID 1 (2 Total Drives)

4GB, Non-ECC, 1333 MHz DDR3

Intel i3 Dual Core Processor 2120 (3.3GHz)

16x DVD +/- RW SATA Optical Drive

Integrated Intel Graphics 2000

Monitor: 2MP, 22" LCD, Resolution:1920x1080 Ratio: 1000:1

Option: External NAS Backup Device 2x1TB (1TB RAID 1)  
a Hardware Component for Datalink.

Note: System Does Not include CAT5 between  
Acquisition Station to the Network Connection.

## CR-X5 Acquisition and Viewing Software

One (1) LAN - Local Viewer / Two (2) Remote Web Viewers

DICOM Send - Multi Destination, DICOM Print

Image Archive

CD/DVD Burning with Mini-Viewer

DICOM Viewer with Robust Imaging Tool Sets

Robust and Advanced Specialty Measurement Tools,

Window Window Level, Rotate, Repeat / Reject

Archive, Image Preview, Magnify, Zoom, Pan, Annotation Edit,

Magnify, Zoom, Pan

L/R Marker, Un-sharp Mask

## Datalink - Disaster Recovery

CR-X5 Systems are provided with a One (1) Year 75GB Disaster Recovery Data Archive Service Plan with an Automatic Backup to a Fully Monitored and Managed HIPPA Compliant Archival and Disaster Recovery Service. Customer's may increase the data plan and/or extend this service after Year One (1) for a fee. Note: Activation of Datalink is automatically activated on all systems prior to shipment.

## CR-X5 Warranty / Support Coverage

Standard Warranty Coverage: One (1) Year CR Parts Replacement and Software Warranty. Internet Access is required for support coverage.

## Workstation PC Warranty and Support Coverage

Computer Acquisition Station and Monitor includes a Three (3) Year Warranty Plan with Next Business Day On-site Support from Lenovo. Help Desk Hours Monday-Friday 9:00AM - 8:00PM Eastern Time, excluding Holidays.