

---

# VIVIX-S Portable, Wired

Portable Flat Panel Detector  
For Digital Radiography



**VIVIX-S Portable, Wired** is a Vieworks' flat panel digital radiography cassette system with 14" x 17" coverage area for general radiographic applications using its unique image processing system and proprietary flat panel detector. With the same size as a film or a CR cassette, it is an excellent solution for upgrading conventional X-ray system. User friendly imaging software, VXvue™ is compliant with DICOM 3.0 standard and provides consistent image quality at a dramatically reduced dose and faster image information with optimized algorithm for each different study. It only takes a few simple steps to acquire and transmit images to the DICOM server through Gigabit Ethernet.

VIVIX

