

Codonics Color Film

Produce stunning color films for cardiac, 3D, SPECT, PET, PET fusion and more with Codonics ChromaVista® Color Film




“Color is a bonus for my referring physicians and is an excellent marketing feature for any practice.”

A. Joseph Borelli, M.D., MRI at Belfair - Bluffton, SC

A4 Size

A Size



- Codonics exclusive Medical Color Matching (MCM™) ensures prints match soft-copy monitors
- Specifically designed for medical applications, our color film is FDA approved and fully archival
- Letter size allows it to fit conveniently into a referring physician report or patient's chart for quick reference
- Ideal for tumor boards and educational purposes




Codonics ChromaVista® Color Film / Transparency




Codonics introduced the world's first DICOM color imager and has been pioneering multi-media imagers for over a decade. We are the industry leader in color imagers with over 25,000 installations worldwide.

Codonics exclusive Medical Color Matching (MCM™) is an advanced image processing feature that adjusts printed output colors to accurately match soft-copy monitors. This process corrects differences in hue, saturation, and intensity between the printer and soft-copy review stations. User's can easily determine the best match for their screen using the automatic MCM bracketing feature which is built into every Codonics imager.



Color prints always match softcopy monitors with our exclusive Medical Color Matching (MCM™) feature. Consistency, clarity and saturation with Codonics color dye-diffusion technology far exceeds that of wax or inkjet printers.



ChromaVista Color Film / Transparencies are ideal for meetings, lectures or legal review.

Specifications:

- Compatible with Horizon and NP/EP multi-media imagers
- Continuous tone color rendering
- Dmin < 0.10 OD, Dmax > 1.80 OD (transmissive)
- Thickness ~ 0.2mm


Available Sizes:

Horizon Imagers

- *A (8.5" x 11")
- *A4 (210 mm x 297mm)

NP/EP Imagers

- A (8.5" x 11")
- A4 (210 mm x 297mm)
- XLA (8.5" x 14")
- XLA4 (210 mm x 356 mm)



For more information contact Codonics at
800-444-1198 or **440-243-1198** and ask for
our Media Sales Department.



17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com

* Horizon color film is a works-in-progress. Contact your Codonics representative for model configurations and media availability. ChromaVista Color Film is available for use in most Codonics imagers, some model and size restrictions apply. Specifications are subject to change without notice. Contact your Codonics representative for the latest information. Copyright © 2005 Codonics, Inc.