

# DARKROOM, VIEWING & FILING SYSTEMS

PRODUCT CATALOG

## Illuminators (View Boxes)

### STANDARD

High quality, cost effective illuminators. Includes upper and lower Film Grips. Internal reflectors offer consistent illumination from two 15-watt fluorescent lamps per viewing area. Use as stand alone or wall mounted. Average light intensity is 2000 nits. 6 foot power cord. UL listed.

#### #VBS-10

14"(W) x 3.5"(D) x 21.75"(H)  
Shipping weight: 17 lbs.

#### #VBD-20

28"(W) x 3.5"(D) x 21.75"(H)  
Shipping weight: 34 lbs.

**SINGLE STANDARD**  
#VBS-10 \$350.00



**DOUBLE STANDARD**  
#VBD-20 \$656.00



## Mammography Illuminator

Light source specifically designed to meet ACR Mammography Viewing requirements. Unit offers two (2) side-by-side film viewing areas for 24" x 30" film. Light intensity: 4,200 nits. Upper and lower film grips. White powder coat finish. Includes four (4) 15-watt fluorescent lamps. UL listed with 5-foot power cord.

#### #VBM-24 | \$779.00

Overall size: 20.5"(W) x 3.5"(D) x 16.75"(H)  
Shipping weight: 17 lbs.



## L.E.D. "X-RAY In Use" Sign

Power-saving LED light technology meets newer compliance demands. Hardwired sign is 11.5"(W) x 7.25" (H). Shipping weight of each sign: 4 lbs.

CAT#	Description	Price
SLU-XL	LED: X-RAY IN USE	\$384.00



## Illuminated "InUse" Door Signs

Lighted sign focuses attention on exam room status. Enhances Radiation Safety by suggesting the risk of accidental entry into an active exam room. NOT a plug-in sign; must be hardwired @ 120VAC. Each sign measures 11.5"(W) x 2.5"(D) x 7.75"(H). Shipping weight of each sign: 4 lbs.

CAT#	Description	Price
SLU-XR	X-RAY IN USE	\$178.00
SLU-CT	CT IN USE	\$265.00
SLU-MR	MRI IN USE	\$178.00



## Standard X-Ray Film Jackets

Choose one or two pockets. Non color-coded (see label selection for sticker styles and colors). Size: 14.5" x 18.25"

including 0.75" scored self expansion. Made of 11pt. manila stock. Black ink copy includes blocks for color-code labels and black loan indicator. Offset double thumbcut for easy access to contents. Double pocket version includes 10" x 15.75" report pocket with 0.75" self expansion. Open side (opens along upper edge).


### SINGLE POCKET

CAT#	QTY Per Carton	Weight Per Carton	Price
JJL-SPS	100	22 lbs	\$189.00



### DOUBLE POCKET


CAT#	QTY Per Carton	Weight Per Carton	Price
JJL-DPS	100	25 lbs	\$222.00



## Mammography Jackets

Pink border identifies as separate from master jacket. Size: 11" x 14", including 0.75" scored self expansion. Pocket size: 8.25" x 11". Illustrated back allows for Dx & patient communication. Offset double thumbcut for easy access to contents. Open side (opens along upper edge).

CAT#	Border Color	QTY Per Carton	Weight Per Carton	Price
JJM-PPM	PINK Border	200	30 lbs	\$368.00
JJM-BBM	BLUE Border	200	30 lbs	\$368.00



## Color-Code Labels

Features permanent adhesive and release coated paper backing. Superior flexibility for long life. Other color systems available. One roll equals 1000 stickers.

CAT#	Description	Rolls	Price	CAT#	Description	Rolls	Price
<a href="#">See List</a>	Any Single Roll	1 roll	\$49.00	<a href="#">See List</a>	Any Single Roll	1 roll	\$40.00
JLL-CNS	0-9 set	10 rolls	\$490.00	JLS-CLS	Set of 10 Rolls	10 rolls	\$400.00

Large digit reverse block



(shown actual size: 1 7/8" x 1 7/8")

Large year



(shown actual size: 1 7/8" x 1 7/8")

CAT#	Color	CAT#
JLL-ORD	red	JLS-CRD
JLL-1GY	gray	JLS-CGY
JLL-2BL	blue	JLS-CBL
JLL-3OR	orange	JLS-COR
JLL-4PR	purple	JLS-CPR
JLL-5BK	black	JLS-CBK
JLL-6YL	yellow	JLS-CYL
JLL-7BN	brown	JLS-CBN
JLL-8PK	pink	JLS-CPK
JLL-9GN	green	JLS-CGN
JLL-CNS	(set of 10)	JLS-CCS
--	aqua	JLS-CAQ
--	gold	JLS-CGD
--	wallet	JLS-CNW
--	white	JLS-CWT

Solid color bar label



(shown actual size: 1.75" x 0.5")

Small year



(shown actual size: 1.75" x 0.5")

CAT#	Year	Color	Rolls	Price
JLL-Y19	2019	Green	1 roll	\$45.00
JLL-Y20	2020	Red	1 roll	\$45.00

CAT#	Year	Color	Rolls	Price
JLS-Y19	2019	Green	1 roll	\$35.00
JLS-Y20	2020	Red	1 roll	\$35.00

## Negative Preservers

Available in three sizes to accommodate all films. All include stock copy lines for: file number, name, date, doctor and remarks. Made of 32lb. gold kraft stock. Thumbcut for easy access to films. Open end (opens along short side edge). Large sized preservers available in red and green.

CAT#	Size	QTY Per Carton	Weight Per Carton	Price
JNS-OEK	8.5" x 10.5"	250	25 lbs	\$242.00
JNM-OEK	10.5" x 12.5"	500	25 lbs	\$262.00
JNL-OEK	14.5" x 17.5"	250	25 lbs	\$204.00

CAT#	Size	QTY Per Carton	Weight Per Carton	Price	Color
JNL-OEG	14.5" x 17.5"	500	25 lbs	\$532.00	Green
JNL-OER	14.5" x 17.5"	500	25 lbs	\$532.00	Red



## Cardboard Insert Mailers

Available in two sizes. Made of 32lb. brown kraft stock. Open end, string & button closure. Chipboard insert for rigidity. Printed with bold black stock copy.

CAT#	Size (inches)	QTY Per Carton	Weight Per Carton	Price
JMM-SBC	11" x 13"	50	7 lbs	\$130.00
JML-SBC	15" x 18"	50	14 lbs	\$158.00



## Green Diamond Border Mailers

Available in two sizes. Made of 11pt. manila stock. Open end, string & button closure. Green diamond border with red. "X-Ray Film, please do not bend."

CAT#	Size (inches)	QTY Per Carton	Weight Per Carton	Price
JMM-SBD	11" x 13"	50	6 lbs	\$158.00
JML-SBD	15" x 18"	50	11 lbs	\$162.00



## E-Z X-Ray Storage Boxes

Allows for the efficient storage of x-ray film files. (Sold in case lots.) Strong and durable with reinforced handles. One piece, easy fold assembly. Can be stacked without crushing. Ideal for archiving files, temporary storage, or moving files. Accommodates jackets up to 14" x 17". Upright: lid opens up. Side down: lid opens forward. (printed legend provided both ways.) Sized for OSHA lift load capacity. Each box O.D.: 6.25" x 16" x 19.5". Ships as 30 lbs. per case. 6 cases or less via FedEx/UPS. 7+ cases via motor freight.

CAT # ETX-70 | \$245.00

(15 boxes per case)



## CD/DVD Sleeves

Standard 5" x 5" CD/DVD protection sleeves. Made of 60lb. white offset paper stock. Window made of clear polyester.

CAT#	QTY Per Carton	Weight	Price
JCD-CSW	1000	14 lbs	\$135.00



## Techno-Tuff™ X-Ray Filing Cabinet

The industry standard cabinets for ARCHIVE or ACTIVE film storage.

Fully assembled, tough, and versatile. All steel construction offers maximum strength. Solid back and side panels for security. 15 3/8" vertical shelf clearance. 12-gauge corner posts for vertical strength. Dividers are 6" apart. Full height dividers add vertical strength and stability. Full depth dividers prevent snagging and ripping of jackets. **Filing cabinets ship via motor freight.**

CAT.#	Height	Depth	Width	Price
608A	84"	18"	24"	\$860.00
618A	84"	18"	36"	\$1015.00
628A	84"	18"	48"	\$1355.00



**#608A | \$860.00**  
Shown in almond



**#618A | \$1015.00**  
Shown in almond



**#628A | \$1355.00**  
Shown in almond

## Receptor Filing Carts

Protect your investment! “CRC” CR-Plate Units have special, impact resistant padding on shelves and vertical surfaces to prevent damage to expensive CR-Plates. Accommodates film jackets, cassettes, and IR plates.

CAT#	Dimensions (W x L x H)	Vert. Clearance	Compartments	Weight	Price
<b>CRC-38</b>	24" x 18" x 42.5"	16.5"	8 - (6" each)	130 lbs	\$1,152.00
<b>CRC-24</b>	24" x 18" x 42.5"	16.5"	4 - (6" each)	100 lbs	\$1,152.00

**Tall Receptor Cart**

**#CRC-38 | \$1152.00**  
Shown in almond







**Mini Receptor Cart**

**#CRC-24 | \$1152.00**  
Shown in almond



**Available Steel Colors**

Platinum  Default	Black  Textured
White 	Almond  Textured

## Standard Wall Bins

Non-shielded steel storage solutions. Will hold receptor plates, film jackets, and caddies. Holds approximately 25 individual file/film jackets. Opening is 22”(W) x 5”(D) x 14”(H).

**#80HP | \$400.00**

HIPAA Version with solid steel no-see front.



**#80AC | \$400.00**

Standard Version with polycarbonate see-thru front.







Corporate Address:

1721 Stewart St. Santa Monica, CA 90404



Reach one of our sales associates at:

1-800-859-6162



For quotes and questions can be directed to:

[info@spectrumxray.com](mailto:info@spectrumxray.com)



Visit us on the web:

[specmedimaging.com](http://specmedimaging.com)

[ppe.spectrumxray.com](http://ppe.spectrumxray.com)