




THYROID COLLARS

PRODUCT CATALOG

Thyro-Guards



Standard

Binding all around means affordable comfort and flexibility.



No Binding

No binding means extra comfort, flexibility and lighter weight.



Mid-Binding

No binding at the top means extra comfort under the chin.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TYR	Regular	\$44.00
TYL	Light Wt.	\$54.00
TYG	Lead-Free	\$64.00

CAT. #	Lead Type	Price
TNR	Regular	\$74.00
TNL	Light Wt.	\$84.00
TNG	Lead-Free	\$94.00

CAT. #	Lead Type	Price
TVR	Regular	\$59.00
TVL	Light Wt.	\$69.00
TVG	Lead-Free	\$79.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TYSR	Regular	\$34.00
TYSL	Light Wt.	\$39.00
TYSG	Lead-Free	\$48.00

CAT. #	Lead Type	Price
TNSR	Regular	\$53.00
TNSL	Light Wt.	\$57.00
TNSG	Lead-Free	\$66.00

CAT. #	Lead Type	Price
TVSR	Regular	\$44.00
TVSL	Light Wt.	\$49.00
TVSG	Lead-Free	\$58.00

Attached Standard

Attached No Binding

Attached Mid-Binding

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TAR	Regular	\$54.00
TAL	Light Wt.	\$64.00
TAG	Lead-Free	\$74.00

CAT. #	Lead Type	Price
TCR	Regular	\$84.00
TCL	Light Wt.	\$94.00
TCG	Lead-Free	\$104.00

CAT. #	Lead Type	Price
TWR	Regular	\$69.00
TWL	Light Wt.	\$79.00
TWG	Lead-Free	\$89.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TASR	Regular	\$46.00
TASL	Light Wt.	\$51.00
TASG	Lead-Free	\$59.00

CAT. #	Lead Type	Price
TCSR	Regular	\$63.00
TCSL	Light Wt.	\$67.00
TCSG	Lead-Free	\$76.00

CAT. #	Lead Type	Price
TWSR	Regular	\$54.00
TWSL	Light Wt.	\$59.00
TWSG	Lead-Free	\$68.00

Attached Thyroid Collars Info

All Attached Thyroid Collars include:

- 1) A snap ring for easy collar attachment.
- 2) Attach anytime with tethered strap.
- 3) Unique 2-in-1 pocket design.

[store collar in one pocket and use other pocket for accessories]

Thyroid collars are 6" tall x 6 1/2" wide, with a 24" strap.



Attached Thyroid option and our 2-in-1 pocket design!

Buckle Closure Add-On

If you prefer buckle closure to regular hook & loop, add "B" to any of the collar numbers above and add \$5.00 to the price.

(Ex.: TYR-B would be \$49.00)



Buckle closure add-on available!

Specialty Thyro-Guards



Pediatric

Binding all around means affordable comfort and flexibility.



Mammography

No binding means extra comfort, flexibility and lighter weight.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TPR	Regular	\$55.00
TPL	Light Wt.	\$60.00
TPG	Lead-Free	\$73.00

CAT. #	Lead Type	Price
TMR	Regular	\$69.00
TML	Light Wt.	\$79.00
TMG	Lead-Free	\$89.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TPSR	Regular	\$53.00
TPSL	Light Wt.	\$57.00
TPSG	Lead-Free	\$66.00

CAT. #	Lead Type	Price
TMSR	Regular	\$53.00
TMSL	Light Wt.	\$57.00
TMSG	Lead-Free	\$66.00

Attached Thyroid Collars Info

All Attached Thyroid Collars include:

- 1) A snap ring for easy collar attachment.
- 2) Attach anytime with tethered strap.
- 3) Unique 2-in-1 pocket design.

[store collar in one pocket and use other pocket for accessories]

Thyroid collars are 6" tall x 6 1/2" wide, with a 24" strap.



Attached Thyroid option and our 2-in-1 pocket design!

Buckle Closure Add-On

If you prefer buckle closure to regular hook & loop, add "-B" to any of the collar numbers above and add \$5.00 to the price.

(Ex.: TYR-B would be \$49.00)



Buckle closure add-on available!


Protective Apparel Fabric Options

NYLON FABRICS

FEATHER RIPSTOP FABRIC *Light Weight Nylon, Filament Reinforced Fabric*



DIAMOND RIPSTOP FABRIC *Medium Weight Nylon, Filament Reinforced Fabric*




NYLON SOLIDS *Liquid Repellent, Light & Supple*



SPECIALTY FABRICS

COLORFUL PRINTS Liquid Repellent, Medium Weight



Prints are shown
SMALLER than
actual size!

TEXTURED FABRIC Textured High Tech Fabrics



Premium Fabric Add-On


NEW Coolwear Fabric

Designed for optimum cooling comfort!

Our new Coolwear fabric add-on provides premium cooling comfort all day long! **This premium fabric can now be added to ANY garment!** The polyester mesh lining features air vented space and stretch-in-width for extreme comfort and fit. This innovative add-on is designed for ultimate performance and comfort.

CP25 - \$25 add-on for all thyroid collars

Features premium comfort
with all day cooling!



Protective Apparel Add-Ons

FREE Embroidery Sewn-In

Select products come with the first line (up to 24 characters) of embroidery **FREE!** Add flexibility with an embroidered patch. Choose from a style below. Add flexibility with an embroidery patch! **Embroidered items are non-returnable & non-refundable.**

Available Styles

JOHN SMITH, MD

John Smith, MD

John Smith, MD

All Caps (default)

Block

Script

#EML2 | \$20 for **second embroidery line**, up to 24 letters.

Thyroid Embroidery Example

Embroidery Colors

- White* (default)
- Black (default)
- Lt. Gray
- Charcoal
- Royal Blue
- Red
- Yellow
- Orange
- Tan
- Brown

Binding Colors

- Black (default)
- Light Gray (default)
- Navy
- Royal Blue
- Red
- Yellow
- Orange
- Tan
- Pink
- Purple
- Green
- Burgundy

Protective Media Guide


MEDIA	BENEFITS	PROTECTION	ATTENUATION	WEIGHT
Regular Lead	Most economical and recyclable	100% Lead	0.5mm Pb Equiv. Protection avg. 97% attenuation over 80, 100, 120 KVP	9 - lbs. (based on medium front protection apron)
Lightweight Lead	Lighter than regular lead and recyclable	Lead-Tin	0.5mm Pb Equiv. Protection avg. 96.8% attenuation over 80, 100, 120 KVP	8 - lbs. (based on medium front protection apron)
Lead Free	Lightest Weight, Lead-Free, Environmentally Friendly	Tungsten-Antimony	0.5mm Pb Equiv. Protection avg. 96.6% attenuation over 80, 100, 120 KVP	7 - lbs. (based on medium front protection apron)



FREE Protective Apparel Tracking and Inventory System!

Exclusive from Techno-Aide!

All new Techno-Aide garments include a scanable QR code for automatic tracking. Upload garment data into your unique ApronHub portal for permanent inventory management. Receive alerts for inspection dates through a customizable format. Retrofit QR code tags onto your existing Techno-Aide or non-Techno-Aide garments. Track inspection history, location, user, etc., through a customizable reports menu. 2D Scanners are available through Techno-Aide or you can use your own 2D-scanner. ApronHub has everything you need to quickly and easily meet the highest QA standards. For more information visit www.apronhub.com today!





Accessories

Tools for Scanning and Labeling Protective Apparel

Stay one step ahead of the inspectors and enhance your QA!
All new garments manufactured by Techno-Aide include a scannable QR

Code Serial Number label for automated or manual tracking. However, we realize that you may have older, unlabeled garments or garments from another manufacturer. To update all your protective garments, we have designed this retrofit labeling kit so that your garments can be tracked utilizing the most up-to-date tracking technology available.

CAT #	Description	Price
AHSCB	2D Scanner with Cable	call
AHSXW	2D Scanner Wireless	call
AHETL	Eyelet Tool	call
AHES50	Silver Eyelets (50)	\$25.00
AHQRC	QRCodeRetrofitLabels(100)	\$68.00
AHDT19	Date Tags: 2019 (50)	\$14.00
AHDT20	Date Tags: 2020 (50)	\$14.00


#AHETL
call for pricing

#AHQRC
\$68.00

#AHSCB/AHSXW
call for pricing

#AHES50
\$25.00

#AHDT19/AHDT20
\$14.00





Corporate Address:

1721 Stewart St. Santa Monica, CA 90404



Reach one of our sales associates at:

1-800-859-6162



For quotes and questions can be directed to:

info@spectrumxray.com



Visit us on the web:

specmedimaging.com

ppe.spectrumxray.com

